# LabgnosticEQA

## **For Laboratories**

Put an end to manually transcribing and typing up External Quality Assurance (EQA) and Proficiency Testing (PT) results, and return schemes directly from your Lab System (LIMS) into your providers' portal.

EQA and PT testing is crucial to ensure your laboratory's testing practices are consistently meeting the highest standard. LabgnosticEQA digitally connects your lab directly with EQA/PT providers, eliminating the current unnecessary manual processes, reducing failures and enhancing efficiency, accuracy, and staff satisfaction.

### **Your Challenge**



Key challenges for labs carrying out EQA/PT reporting:

- Time and cost inefficiency
  An enormous amount of time is spent on manual
  EQA/PT processes which creates unnecessary time
  and cost inefficiencies.
- Reduces lab capacity for patient testing
  Valuable staff time spent on EQA/PT reporting could be redeployed to turning around patient samples and other business-critical laboratory tasks.
- Frustrating failed returns and low EQA/PT scores
  Filling out forms to onboard samples to analysers and
  manually transcribing results leads to delays, errors and
  inaccurate EQA/PT reports, contributing to lower EQA/PT
  scores or failures.
- Lower staff satisfaction

  The administrative burden of EQA/PT reporting is labourintensive and highly pressurised which can lead to
  dissatisfaction among laboratory staff who are
  responsible for maintaining the high EQA/PT scores.








Key advantages of using Labgnostic EQA for EQA/PT reporting:

- Removes manual processing
  LabgnosticEQA automates and streamlines the manual process of booking-in, onboarding, and resulting samples. It also eliminates the need to check for transcription errors before returning results.
- Redeploys time to the lab
  The digital workflow reduces the time laboratory staff spend on EQA/PT reporting, increasing operational efficiency by giving time back to patient testing and other business-critical activities.
- 3 EQA/PT reporting mirrors patient sample handling
  The seamless integration of your lab system (LIMS) with its
  EQA/PT providers enables staff to handle EQA/PT samples
  like normal patient samples. The process is more familiar and
  ensures that the lab's results accurately reflect the lab's
  quality.
- Efficient onboarding
  Our dedicated deployment team ensures a smooth transition for new labs, completing the setup in just a few weeks.



# How does LabgnosticEQA work?

Receive sample requests and send results instantly via a direct digital connection, with the added benefits of specimen tracking and powerful customization features.





"Using Labgnostic for our EQA/PT schemes means we now have zero requesting errors, transcription errors, and failed returns for lateness. The number of staff needed to perform EQA has been reduced and staff have actually said that they feel less stressed by the new process and retaining our great EQA scores."

Adam Tonge, IT Coordinator in the Department of Biochemistry at Manchester Royal Infirmary

## **Benefits**



#### **Reflective EQA Scores**

Choose LabgnosticEQA for the certainty that your EQA/PT results reflect the quality of your lab's testing and reporting processes.



#### **Happier Labs**

Ensure staff feel more satisfied that the EQA/PT scores reflect the quality of their work and give them the capacity to focus on what matters.



#### **Time and Cost Saving**

Increase operational efficiency by cutting the time and cost associated with EQA/PT reporting and optimise valuable laboratory resources for patient testing.



#### **No Unnecessary Failures**

Improve compliance by eliminating transcription errors and failures for late returns.

### **Our Partners**



**RCPAQAP Australia** 



**UK NEQAS Birmingham Quality** 



**UK NEQAS Haematology** 



UK NEQAS Immunology, Immunochemistry & Allergy





With 95% penetration in the UK NHS market, X-Lab processed more than 500 million tests for UK citizens during Covid. Labgnostic, our lab-to-lab interoperability solution, is the world's first system-agnostic laboratory exchange network. And LabgnosticEQA delivers services for EQA/PT providers worldwide. X-Lab operates throughout the UK, Europe, Australia, and the USA.



www.x-labsystems.com



info@x-labsystems.co.uk



@labgnostic

3000



X-Lab Level 7, 330 Collins Street, Melbourne, Vic,



time spent on manual EQA/PT processing

reduction in staff

Labs have reported a