

Labgnostic PT

For Laboratories

Put an end to manually transcribing and typing up Proficiency Testing (PT) results, and return schemes directly from your Lab Information System (LIS) into your providers' portal.

PT testing is crucial to ensure your laboratory's testing practices are consistently meeting the highest standard. LabgnosticPT digitally connects your lab directly with PT providers, eliminating the current unnecessary manual processes, reducing compliance failures and enhancing efficiency, accuracy, and staff satisfaction.

How does LabgnosticPT work?

Receive sample orders and send results instantly via a direct digital connection, with the added benefits of specimen tracking and powerful customization features.



"Using Labgnostic for our PT schemes means we now have zero requesting errors, transcription errors, and failed returns for lateness. The number of staff needed to perform PT has been reduced and staff have actually said that they feel less stressed by the new process and retaining our great PT scores."

Adam Tonge, IT Coordinator, MRI Hospitals



Benefits



Zero errors

Using Labgnostic to submit PT returns will eliminate ordering and results errors, as well as failed returns.



Saves time

Digitally receiving and submitting results reduces staff time and stress. Staff can be redeployed to meet growing workload demands within the lab.



No late returns

Labgnostic minimises transcription errors. It reduces pre- and post-analytical errors and eliminates the possibility of late returns.

100%
of labs

Have reported a reduction in staff time spent on manual PT processing

Your challenge



Time and cost inefficiency

An enormous amount of time is spent on manual PT processes which creates unnecessary time and cost inefficiencies.




Reduces lab capacity for patient testing

Valuable staff time spent on PT reporting could be redeployed to turning around patient samples and other business-critical laboratory tasks.



Frustrating failed PT returns and poor CLIA compliance


Filling out forms to onboard samples to analysers and manually transcribing results leads to delays, errors and inaccurate PT reports, contributing to lower PT scores or failures.



Lower staff satisfaction

The administrative burden of PT reporting is labor-intensive and highly pressurised which can lead to dissatisfaction among laboratory staff who are responsible for maintaining the high PT scores.

Our solution



Removes manual processing

Labgnostic automates and streamlines the manual process of booking-in, onboarding, and resulting samples. It also eliminates the need to check for transcription errors before returning results.



Redeploys time to the lab

The digital workflow reduces the time laboratory staff spend on PT reporting, increasing operational efficiency by giving time back to patient testing and other business-critical activities.



PT reporting mirrors patient sample handling

The seamless integration of your lab information system (LIS) with its PT providers enables staff to handle PT samples like normal patient samples. The process is more familiar and ensures that the lab's results accurately reflect the lab's quality.








Efficient onboarding

Our dedicated deployment team ensures a smooth transition for new labs, completing the setup in just a few weeks.

About X-Lab

With operations in the UK, Ireland, EU, USA, New Zealand, and Australia, X-Lab is a global company. Its products have saturated the UK market and processed more than 500 million tests during the pandemic. X-Lab services the lab2lab interoperability market, the remote order comms market and facilitates proficiency testing for PT providers and labs worldwide.

Contact Us

-  www.labgnostic.com
-  hello@labgnostic.com
-  X-Lab Global
-  @Labgnostic
-  Los Angeles, CA